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Now I would like to talk about the actual handling of radiation generating devices.

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This time, we will explain what a radiation generator is, the role of beamlines in using radiation generators, the role of beamlines in safety measures, and finally the procedures for using a radiation generator.

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A radiation generator is generally defined as any device that generates radiation, and is generally interpreted to include ordinary X-ray equipment; however, the Act on the Regulation of Radioactive Isotopes, etc. defines it as follows, and does not include ordinary X-ray equipment.

The following types of devices are defined as radiation generating devices: devices that generate electron beams or X-rays (limited to those with an energy of 1 MeV or more) and other radiation by accelerating charged particles such as electrons and protons using electric or magnetic fields, and with a maximum dose equivalent rate exceeding 600 nSv/h at a position 10 cm away from the surface of the device.

① Cyclotron ② Synchrotron ③ Synchrocyclotron ④ Linear accelerator ⑤ Betatron ⑥ Van de Graaff accelerator ⑦ Cockcroft-Walton accelerator ⑧ Transformer accelerator ⑨ Microtron ⑩ Plasma generator (device capable of achieving the critical plasma conditions in DT nuclear reactions, limited to those dedicated to DT nuclear reactions)

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
Following the Physics Prize, the 2024 Nobel Prize in Chemistry was also awarded for an achievement related to AI technology.

Radiation generators are widely used in the development of semiconductors, protein structural analysis, and drug development, and we believe that databases that consolidate various analytical information accumulated to date have been utilized.

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
According to data from the 2023 edition of the Nuclear White Paper, a total of 1,747 radiation generating devices have been approved for use in Japan, of which approximately 75% are owned by medical institutions and are linear accelerators (linacs). Others are installed by research institutes such as the RIKEN, educational institutions such as Tohoku University, and private companies.

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
The acceleration energy is 8 GeV, which is greater than that of the APS in the US and the ESRF in Europe.

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In Japan, there are RIKEN's Spring-8 8GeV in Harima Science Garden City, Hyogo Prefecture, NewSUBARU 1.5GeV in Hyogo Prefecture, Hiroshima University's Synchrotron Radiation Research Institute Hisor 700MeV, Saga Prefectural Kyushu Synchrotron Light Research Center SAGA-LS 1.4GeV, Ritsumeikan University Research Organization for Science and Technology SR SENTER Rits-SR 575MeV, Aichi Synchrotron Light Center 1.2GeV, Institute for Molecular Science Extreme Ultraviolet Research Facility UVSOR 750MeV, High Energy Accelerator Research Organization Synchrotron Radiation Experimental Facility Photon Factory PF Ring 2.5GeV, Advanced Ring PF-AR 6.5GeV, and other research, educational, and medical facilities that use accelerators scattered from Hokkaido to Kyushu. Select the facility

you want to use depending on your purpose.

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In addition, the NanoTerrace, a new 3 GeV high-brilliance synchrotron radiation facility, has begun operation in Sendai, Miyagi Prefecture.

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The Ministry of Education, Culture, Sports, Science and Technology website has videos introducing radiation generating equipment facilities. Due to time constraints, I will only introduce one from KEK. If you are interested, please watch it later.

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To get an idea of what goes on inside a synchrotron radiation facility, check out this video about RIKEN's X-ray Free Laser.

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
This is an aerial photograph of Spring-8. Comparing it with the following image will give you an idea of the facilities and functions inside the building.

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
Electrons fired from the electron gun installed at SAKURA in the lower right are accelerated to 8 GeV by SACLA's linear accelerator and injected into a storage ring with a circumference of approximately 1,500 m, where they are rotated while maintaining their energy, generating synchrotron radiation using bending magnets and an insertion light source. The generated synchrotron radiation is guided through beamlines to experimental hatches installed inside and outside the storage ring, and is used for a variety of experiments.

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
Synchrotron radiation is a high-intensity white electromagnetic wave that is generated in the tangential direction by bending the orbit of electrons accelerated to nearly the speed of light using an electromagnet.

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Synchrotron radiation is characterized by its brightness, energy, and directionality. This is from Spring-8, which is 1 billion times brighter than regular X-ray equipment, and has energies greater than infrared, ranging from soft X-rays of 170 eV to hard X-rays of 300 keV, as well as gamma rays up to 2.9 GeV.

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
Next, I would like to discuss the role of beamlines in utilizing radiation generators.

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
In order to extract the type and energy of radiation appropriate for the user's experimental purpose, the light is focused, shaped, and monochromatized using spectrometers and mirrors installed in the beamline.

The X-rays that you will be handling are produced by radiation generating devices such as linear accelerators and synchrotrons. Therefore, since you will be using radiation generating devices, you are eligible to take this education and training.

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Here, we will introduce the beamlines at Spring-8. As you can see, there are 62 beamlines available, and they are used for experiments in a wide range of fields, including soft X-ray solid spectroscopy, soft X-ray photochemistry, white X-ray diffraction, photoelectron spectroscopy/microCT, medicine/imaging, and crystal structure analysis.

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This is the experiment station. The radiation is introduced into the experiment hatch in the experiment hall, where it is irradiated onto samples and measurements are taken.

This photo shows the scene where the samples and measurement equipment are being set up. After setting up is complete, we make sure that no one is in the experimental hatch, close the entrance door, turn on the interlock, and start irradiation and measurement.

If the emergency stop button is pressed, all beamlines may enter an emergency shutdown. It is important to make sure that no samples or instruments are set up and no people are present in the experiment hatch.

There is an emergency stop button near the entrance to the experimental hatch.

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
In cases where ultraviolet light or radiation, which have little effect on the human body, is used, no experimental hutch is provided. For example, there is the Synchrotron Radiation Facility (UVSOR) at the Institute for Molecular Science.

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This is an overview of the Kyushu Synchrotron Light Research Center. The maximum energy is 1.4 GeV, which is smaller than the 8 GeV of Spring-8, but it is used for experiments in a variety of unique fields.

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
This is the experimental hall of the Kyushu Synchrotron. Experimental hutches are set up according to the type and energy of radiation used.

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X-rays generated from synchrotron light are used in a variety of experiments, taking advantage of the interaction between the X-rays and materials. By measuring the scattered rays, crystal structures can be analyzed, and local structures and chemical bonding states can be determined by imaging based on the penetration conditions, while elemental analysis can be performed by measuring fluorescent X-rays.

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
The characteristics of radiation obtained from synchrotron radiation include high brightness, high directionality, and variable energy. It enables trace element measurement that could not be performed with normal X-ray devices, shortens measurement time, and enables high-resolution measurement.

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Synchrotron radiation contributes to the advancement of various analyses that utilize the interaction between X-rays and materials.

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Here are four of the recent results obtained using synchrotron radiation. This is the development of a protein crystal structure analysis method.

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This is the discovery of the three-dimensional arrangement of dopamine atoms in diamond by photoelectron holography.

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This is the successful determination of the ordered growth and structure of nano-composite crystalline thin films by ultra-thin film structure X-ray diffraction. It is expected that in the future, this technology will enable the storage and transportation of liquefied gas to be changed from the conventional heavy cylinders to lighter, easier-to-handle boxes.

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
This project used micro-CT to elucidate the formation and evolution of the asteroid Ryugu.

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This is a heavy ion cancer treatment device at the National Institutes for Quantum and Radiological Science and Technology. By using heavy carbon beams obtained from synchrotron radiation, it is possible to maximize the therapeutic effect on the target tissue while reducing the side effects that occur when radiation is absorbed by tissues other than the affected area, as occurs with conventional X-rays and gamma rays.

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
Next, I would like to talk about the role of beamlines in safety measures.

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

To prevent users from being exposed to radiation and equipment damage, and to maintain safety by maintaining a vacuum in the storage (light source) ring, beamlines are equipped with radiation shielding measures and interlock systems.

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
The synchrotron radiation extracted from the storage ring contains high-energy gamma rays and neutrons, which are blocked by absorbers and shielding walls.

In addition, a gamma ray stopper is installed after the spectrometer to adjust the energy of the radiation introduced into the experimental hut.

Furthermore, an interlock system consisting of a radiation shutter and a downstream shutter is installed in case of an emergency inside the experimental hut.

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
Synchrotron radiation facilities are required to designate and manage as controlled areas the radiation dose at the boundary of their workplaces where it is kept below $250\ \mu\text{Sv}$ per three month period, and any area where the radiation dose exceeds 1.3mSv per three month period.

Those who wish to use the facility must be registered as radiation workers at Kumamoto university and must have undergone education and training and a health check.

When using the synchrotron radiation facility, please bring your glass badge issued by Kumamoto University.

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
Users must follow the instructions of the radiation protection supervisor and safety manager at each synchrotron radiation facility, including time restrictions on access to restricted areas and experimental locations. The radiation dose limit at the usage location is set at 1 mSv per week.

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Neutrons are produced by nuclear reactions (γ -n) between electromagnetic waves of 6 MeV or higher and matter.

These neutrons activate the concrete and metal of the shielding walls and the materials of the beamline.

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Finally, we will explain the procedures for using synchrotron radiation.

放射線作業登録申請書(2022年度)
Radiation Worker Application Form
(For External User / Supplier)

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This is the application form for Spring-8. Recently, applications have been accepted online. Even though it is an online application, you cannot apply on your own.

Before registering online, please confirm the details you have entered with the relevant administrative department of your department or department.

諸規定を遵守し、放射線作業に従事することを誓います。
I will take the classes at the request of the laws and the RIKEN rules, and understand the contents of them. I hereby agree to do the radiation work under the condition provided by them, and promise to keep them.

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所属機関承諾・証明欄 Approved and authorized by affiliated organization as below:

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所属機関名称
Name of affiliated organization/company

住所
Address

代表者職・氏名
Name and title of the representative for the organization above

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The details will be checked by the head of the department in charge and the radiation protection officer. If there are no discrepancies with the facts, approval will be processed via the web system.

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Please strictly observe these precautions.

- 1 You must be registered as a radiation worker permitted by your affiliation in order to engage with radiation work at this campus.
- 2 If a radiation protection supervisor of this campus requested to submit your records of exposure dose, medical examinations, and radiation training, please submit these promptly.
- 3 When you engage with radiation work at this campus, please bring a dosimeter belonging to your affiliation.
- 4 Before you start radiation work, please take safety trainings (e-learning or DVD training) and understand radiation hazard control regulations of this campus.
- 5 Please make efforts for security for radiation related works.

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Conditions for using synchrotron radiation.

1. The individual must be a "radiation worker" employed by Kumamoto university.

2. 2. Obtain approval from the head of the department or division to which you belong and the radiation protection officer.
3. Have undergone a health check before entering the controlled area and regularly (January and July).
4. Have taken the required education and training.
5. If you plan to use synchrotron radiation, you will need to bring along a "Neutron Wide Range (NS)" instrument. Therefore, if you plan to use synchrotron radiation, contact the department that will issue the measuring device immediately.

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I used these documents as reference.
Thank you for taking this course.